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1 x CR2032 battery The device can emit effective output values above 10 mA, ave aged over every five-second interval ፟
大 Application part, type BF Dispose of packaging in an environmentally friendly manner Certification symbol for products that are exported to the Russia Federation and CIS countries

Please read these instructions carefully before using the device.

LOT number

Expiry date

Manufacturer

. Getting to know your device 1.1 What is the EMS-pad?

he innovative EMS Sixpack Pad is a self-adhesive EMS device for muscle training and axation. The device can be applied with accuracy thanks to the use of EMS technology Made from medical-quality silicone rubber, the elegantly designed EMS pad is extremely slim and flexible and adapts perfectly to the area to be treated. The high-tech circuit pro vides energy-efficient treatments of 20 minutes each (up to 20 treatments per battery 1.2 What is EMS and what can it do? Electrical muscle stimulation (EMS) is a widespread and generally recognised nethod and has been used in sports medicine and rehabilitation for years. In sports and

fiftness, EMS is used to complement conventional muscle training, to increase the per-formance of muscle groups and to adjust physical proportions to achieve the desired aesthetic results. There are two different types of EMS application: One is for targeted strengthening of the muscles (activating application), and the other is to achieve a elaxing, restful effect (relaxing application). The activating application involves: Muscle training to increase endurance and/or Muscle training to support the strengthening of specific muscles or muscle groups, and

to achieve the desired changes to physical proportions

The relaxing application involves: n for easing muscle tension ration of muscle regeneration after high muscle performance (e.g. after a

EMS devices work by passing electrical currents over the skin. The gel film is used as a transfer medium and is subject to natural wear and tear. The gel film must be replaced if it stops providing sufficient contact, as this will prevent the EMS pad from sticking to the skin. If it is not replaced, the partially increased current density could irritate the skin 1.3 Included in the delivery 1 x EMS-pad x self-adhesive gel film set 1 x 3V CR2032 battery Nehmen Sie die Batterie aus dem Artikel heraus, wenn diese erschönft ist oder Sie

 1 x self-adhesive gel film set and 1 x battery, item no.: 647.16 When necessary, please contact Customer Services at the address provided in good time N IMPORTANT! Beschädigung

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· On people.

• For the intended purpose and as specified in these instructions for use. Any improper • For hygienic reasons, the device may only be used on one person. • If the device does not work properly, or you feel unwell or experience pain, stop us-Use only with the original accessories supplied, which can be re-ordered. Failure to do

 Switch off the EMS pad before attaching or removing it. Do not adapt the EMS pad in any way (e.g. by cutting or trimming it). This will produce The device is not a substitute for medical consultation and treatment Consult your doctor first if you are experiencing any pain or are suf a higher current density and can be dangerous.Do not use the device whilst using other devices that transmit electrical impulses into your body.

 Do not use whilst asleep, driving a vehicle or operating machinery Do not use whilst undertaking any activity where an unexpected reaction (e.g. strong muscle contractions even at low intensity) could be dangerous.
 Ensure that no metallic objects (e.g. belt buckles or necklaces) come into contact with the electrodes during stimulation. If you are wearing jewellery or have piercings in the area to be treated (e.g. a navel piercing), these must be removed before using the derice. Failure to do so could result in spot burns.

 Keep the device away from children.
 Pre-treating the skin with greasy creams or lotions is not recommended. Do not use adhesive tapes, bandages or any aid other than the designated gel film to apply the EMS pad to your skin. Otherwise the unequal distribution of the gel film could lead to skin lesions. The device should be used under medical supervision where possible. · Stimulation on the same area for a long period of time could irritate the skin. Please

nsult a dermatologist if skin irritation persists. ring the initial few minutes, use the device while sitting or lying down to minimise the risk of accidental injuries as a consequence of isolated cases of vagal responses (feeling of faintness). If you feel faint, switch off the device immediately, lie down and ort the legs in an elevated position (approx. 5 – 10 min) Notes on handling batteries

If your skin or eyes come into contact with battery fluid, flush out the affected areas with water and seek medical assistance. • 🛧 Choking hazard! Small children may swallow and choke on batteries. Store the Observe the plus (+) and minus (-) polarity signs. . If a battery has leaked, put on protective gloves and clean the battery compartment

with a dry cloth. Protect the batteries from excessive heat. A Risk of explosion! Never throw batteries into a fire. Do not charge or short-circuit batteries. If the device is not to be used for a long period, take the batteries out of the battery

identical or equivalent battery types only Do not use rechargeable batteries. · Do not disassemble, split or crush the batteries

 This unit is not intended for use by children or persons with physical, sensory (e.g. in-sensitivity to pain) or mental impairments or those who do not know how to use it due to a lack of experience or knowledge, unless they are supervised by a person respon-sible for their safety or have received instructions from them on how to use the unit. • Remove the battery from the item if it is empty or you do not intend to use the item for a prolonged period of time. This prevents damage which could occur due to leakage Before using the device, consult your doctor if any of the following uired, clean the battery and device contacts before inserting. Risk of overh • Serious illnesses, in particular if you suspect or have been diagnosed with high blood If the EMS pad is damaged or you have any doubts, do not use it and contact your re-

pressure, a blood coagulation disorder, propensity to thrombo-embolic conditions or tailer or Customer Services at the address provided. Switch the device off immediately if it is faulty or not working properly Do not attempt to open and/or repair the EMS pad yourself. Repairs may only be carried out by Customer Services or authorised suppliers. Failure to comply will result in Any sensory impairment that reduces the feeling of pain (e.g. metabolic disorders). voiding of the warranty. The manufacturer is not liable for damage resulting from improper or careless use.

 If you suffer from persistently irritated skin due to long-term stimulation at the same. Remove the batteries from the device before cleaning. Clean the device after use with a soft, slightly damp cloth. If it is very dirty, you can also Total the cloth with a mild soapy solution.
 Ensure that no water gets into the device. If this happens, only use the device again

once it has fully dried out. 1 m) to shortwave or microwave devices (e.g. mobile phones), as doing so can result Do not use any chemical or abrasive cleaning agents

• Do not expose the device to direct sunlight or high temperatures.
• Protect the device from dust, dirt and humidity. Never immerse the device in • Re-apply the protective film to the gel film after use.

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ENGLISH

avoid damage to health, we strongly advise against using the EMS

. If you have a known or acute cardiac arrhythmia, or disorders of the heart's impulse

After an operation, if strong muscle contractions could affect the healing process.
 If you have had a back operation in the past and the back pain is due to the surgery

On acutely or chronically diseased (injured or irritated) skin (e.g. inflamed skin - wheth

In humid environments (e.g. in the bathroom) or when bathing or showering.

er painful or not, reddened skin, rashes, e.g. allergies, burns, bruises, swellings, both open and healing wounds, and post-operative scars where the healing process could

The device must never be used near to the heart. The stimulation elec-

trodes must not be placed on any part of the front ribcage (where the ribs

and breastbone are located), especially not on the two large pectorals. This can increase the risk of ventricular fibrillation and induce cardiac arrest.

use can be dangerous.

2. Important notes

Instructions for application

If you have a high temperature (e.g. > 39 °C).

or an accident.

If you suffer from metabolic arthritis.

Do not use after consuming alcohol.

ecurrent malignant growths.

If you are receiving medical treatment.

PRECAUTIONS:

Per uso esterno.

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If you suffer from a seizure disorder (e.g. epilepsy)

. for acute or chronic diseases of the gastrointestinal tract.

In the neck area/carotid artery area/above the neck

Do not use if connected to a high-frequency surgical device.

ained chronic pain in any part of the body

If you are already having stimulation treatment for existing condition

. Hold the device away from sources of heat and do not use it in close proximity (approxi-

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If you have implanted electrical devices (e.g. a pacemaker).

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• If you do not plan to use the device for a long period of time, we recommend that you

npartment lid anti-clockwise using a co

Use a pointed object to loosen the lid.
 Place the battery in the battery compartment so that the positive terminal (+) fac-

. Remove one of the protective films carefully.

film does not protrude over the EMS pad.

• Carefully and slowly remove the protective film, making sure that the two parts of the

adhesive gel film remains on the same side of the protective film. Ensure that the ad-

hesive gel film does not present any damage because if it is damaged or not smooth

The EM 20 Sixpack was specially developed for building and toning abdominal muscles

Stick the EMS pad onto the part of the body to be treated (see image)

1)

2. Carefully apply the gel film onto the EMS pad and then carefully

peel away the protective film. Ensure that the edge of the gel

Gel film that has been applied slightly askew will not affect the

1. Press and hold down the "+" button until you hear a long beep

2. Press the "+" button again to start the Level 1 stimulation. Keep

To turn off the EMS pad, press and hold down the "-" button until you hear a long beep. The device switches off automatically af-

3. Press the "-" button to reduce the intensity.

For optimal results we recommend a 4-week training plan. Please note the following points

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Re-apply the protective film to the gel film after use.

3. Initial use

I Inserting the batter

0

4. Close the lid by turning it clockwise.

3.2 Attaching the gel film

0

2 Increase the intensity level on a daily basis

3. A higher intensity level can lead to better results. However, the stimulation should

Wait at least 30 minutes after the stimulation programme is complete before using the

• The device will only work if there is sufficient contact with the skin. Stimulation

Replace the battery if it is low (insufficient stimulation, no beeping noises).
 For optimum stimulation we recommend that the gel film is replaced when the

In the event of continuous beeping, shut down the device. Clean the gel film and try to restart the program. If the continuous beeping persists, replace the gel film.

If the adhesive capability of the gel film decreases, replace immediately. Do not use the EMS pad again until the gel film has been replaced. Otherwise the unequal

distribution of the gel film could lead to skin lesions.

In the event of unpleasant current peaks, switch off the EMS pad and use it in a dif-

• The "Smart Contact Detection" automatically triggers an audible alarm if 2 electrodes

400 μs 400 μs 400 μs

2.5 s 2.5 s 1.5 / 10 s **

never be uncomfortable.
4. Complete the full 20-minute programme each time you use the device.

nulses will not be transmitted if there is insufficient contact

Do not use the EMS pad on the same area more than once per day.

Switch off the EMS pad before removing it.

ome into contact with each other.

Pulse frequency 2 Hz. 4 Hz

* There are 3 individual rounds in cycle 3.

** The off-time is 10 seconds for every 3rd round in cycle 3.

ers. You are legally required to dispose of the batteries.

Pb = Battery contains lead, Cd = Battery contains cadmium

contact the local authorities responsible for waste disposa

EM 20 EMS Sixpack Pad

reatment time 20 minutes, the device then switches off automatically

Storage conditions -10 °C - 50 °C (14 °F - 122 °F) with a relative humidity of

66 g (without battery), 69 g (with battery)

If the device is not used according to the instructions specified, perfect functioning can

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We reserve the right to make technical changes to improve and develop the product.

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Press and hold down the "+" button agai - Check that the battery is inserted correctly and try t witch on the device again. The battery could be flat. Replace the battery and try t switch on the device again. ally during use. witch off the device. Then press and hold down the button for 3 seconds to switch on the device again. The batteries could be low and must be replaced. You are unable to feel any - Press and hold down the "+" button again for Ensure that the gel film does not protrude from the d vice and try again.

- Gradually increase the intensity level. The skin at the area of - Stop treating this area immediately. Wait until the skin eatment is red. has returned to its normal colour. Please consult a dermatologist if skin irritation persists. The intensity suddenly be- - Replace the battery and try again.

8. Warranty/service

In case of a claim under the warranty please contact your local dealer or the local reprecase of returning the unit please add a copy of your receipt and a short report of

The following warranty terms shall apply: 1.The warranty period for BEURER products is either 3 years or- if loner- the country specific warranty period from date of Purchase.

a case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice. !. Repairs (complete unit or parts of the unit) do not extend the warranty period. 3. The warranty shall not be valid for damages because of

 a. improper treatment, e.g., nonobservance of the user instructions.
 b. repairs or tampering by the customer or unauthorised third parties.
 c. transport from the manufacturer to the consumer or during transport to the ser-For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE and tear (cuff, batteries etc.). (Waste Electrical and Electronic Equipment). If you have any questions, please

d. The warranty shall not be valid for accessories which are subject to normal wear 4. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

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